




THD STERILE SINGLE-USE DIAGNOSTIC ANOSCOPE



Use

The THD Diagnostic Line of sterile diagnostic anoscopes is made up of a wide range of devices available in various diameters and openings to satisfy the various clinical and pathological requirements of anoscopy practices in a surgical and out-patient environment.

General description


The anoscope is constructed of medical-grade, transparent and shockproof thermoplastic material. The anoscope can be used with the THD LED2000 illuminator and THD REVOLUTION (illuminator function) or with an operating theatre illuminator. The dedicated shaped

fibre-optic tip which is available as an accessory allows effective illumination of the operating area as the light is aimed directly at the tip.

Each anoscope complete with its own introducer is individually packed in a sterile double bag

Features

- Made of medical-grade polycarbonate
- Complies with the biocompatibility requirements established by ISO 10993
- Atraumatic insertion profile
- Ready to use
- Class I CE device
- Individually packed in polyethylene bags in a practical "Gravity Feed" box
- Latex-free
- Complete range of diameters, lengths and openings
- Diameter quickly identifiable through the handle colour
- Compatible with the THD LED2000 illuminator and THD REVOLUTION
- Sterile
- Adaptable to other light sources



		Paediatric	Medium	Large
Handle colour		Orange	Blue	Green
d	Internal lumen	12 mm	19 mm	23 mm
D	Outside diameter	15 mm	22 mm	26 mm
L	Total length	65 mm	90 mm	90 mm

THD STERILE single-use diagnostic anoscope		
TRUNCATED		BEVELLED
Code	Qty	Denomination
800011 - 90	90	Sterile Paediatric Truncated Anoscope
800012 - 90	90	Sterile Paediatric Bevelled Anoscope
800014 - 45	45	Sterile Medium Truncated Anoscope
800015 - 45	45	Sterile Medium Bevelled Anoscope
800017 - 45	45	Sterile Large Truncated Anoscope
800018 - 45	45	Sterile Large Bevelled Anoscope

THD360049_EN_Rev00_20180503_THD_Diagnostic Ano Sterili.docx