




THD SINGLE-USE DIAGNOSTIC ANOSCOPE



Use

The THD Diagnostic Line of diagnostic anoscopes is made up of a wide range of devices available in various diameters and openings to satisfy the various clinical and pathological requirements of anoscopy practices in an out-patient environment.


General description

The anoscope is constructed of medical-grade, transparent and shockproof thermoplastic material. The anoscope can be used with the THD LED2000 illuminator and THD REVOLUTION (illuminator function) or with an operating theatre illuminator. The dedicated shaped fibre-optic tip is available as an



accessory which allows effective illumination of the operating area as the light is aimed directly at the tip. Each anoscope complete with its own introducer is individually packed.

Features

- Made of medical-grade polycarbonate
- Complies with the biocompatibility requirements established by ISO 10993
- Atraumatic insertion profile
- Ready to use
- Class I CE device
- Individually packed in polyethylene bags in a practical "Gravity Feed" box
- Latex-free
- Complete range of diameters, lengths and openings
- Diameter quickly identifiable through the handle colour
- Compatible with the THD LED2000 illuminator and THD REVOLUTION
- Adaptable to other light sources



		Paediatric	Medium	Large
Handle colour		Orange	Blue	Green
d	Internal lumen	12 mm	19 mm	23 mm
D	Outside diameter	15 mm	22 mm	26 mm
L	Total length	65 mm	90 mm	90 mm

THD single-use diagnostic anoscope		
TRUNCATED		BEVELLED
		
Code	Qty	Denomination
800020 - 120	120	Paediatric Truncated Anoscope
800021 - 120	120	Paediatric Bevelled Anoscope
800023 - 90	90	Medium Truncated Anoscope
800024 - 90	90	Medium Bevelled Anoscope
800026 - 90	90	Large Truncated Anoscope
800027 - 90	90	Large Bevelled Anoscope

THD360050_EN_Rev00_20180503_THD_Diagnostic Ano NS.docx