




# THD LIGHTSCOPE ANOSCOPE



## Use

The LightScope anoscope is a DISPOSABLE medical device deemed for investigation of the anal canal.

## General Description

The LightScope anoscope can be used with the specific THD optical fiber end tip P/N 880004/880012. The LightScope anoscope has been specifically designed to be used with the THD optical fiber end tip and it has a flat white end to reflect the light and have a clear and detailed vision of the inspection field.

Available with truncated or beveled tip.

## Features

- Made of a medical grade Polycarbonate
- Comply with the biocompatibility requirements of the ISO 10993 standard
- Specifically designed mandrel
- Two ends available (truncated and beveled)
- Ready to use
- Light (40 g/0.09 pound)
- CE Class I
- Individually packaged
- Latex Free



Advantages

- Effective lightning
- Compatible with the THD end tip (P/N. 880004/880012)

Dimensional Features	
THD LIGHTSCOPE ANOSCOPE	
Diameter (B)	25 mm (0.98 in)
Length (A)	70 mm (2.76 in)

THD LightScope Anoscope		
TRUNCATED		BEVELLED
Code	Qty	Denomination
800102 - 10	10	THD LightScope Ano Truncated
800101 - 10	10	THD LightScope Ano Bevelled

THD360130\_EN\_Rev00\_20190424\_THD\_Lightscope\_Ano\_NS.docx