




THD SURGY MINI LIGHTSCOPE



Use

The SURGY LightScope anoscope is a SINGLE-USE medical device indicated for anorectal canal analysis and surgery.

General description

The SURGY LightScope anoscope combines a practical and lightweight anoscope and a powerful integrated LED light source in one single-use device in order to give you a clear and detailed view of the inspection area.

The SURGY LightScope anoscope is ready for use and does not require any other equipment for use.

The light source is activated simply by pulling out a tab.


Features

- Made of medical-grade polycarbonate
- Complies with the biocompatibility requirements established by ISO 10993
- **Self-light, with a permanent integrated light source**
- Atraumatic insertion profile
- Extremely compact
- Ready to use
- Lightweight (40g)
- Class I CE device
- Individually packed
- Latex-free



**Advantages:**

- Does not require an additional light source
- Does not require preparation before use
- Eliminates the risk of cross-contamination
- Compact and easy to handle
- Mini-invasive and atraumatic



	Surgy mini
E	27 mm
LF	122 mm
I	23
L	85



After use, dispose the THD Light-Scope rectoscope like a biological waste according to the local laws.



If the THD Light-Scope rectoscope with LED needs to be disposed before it has been used (i.e. because it is expired or damaged), do not dispose as general waste. It contains electric components.

THD Surgi Mini LightScope			
Code	Qty	Trade Name	Denomination
800108	10	THD LIGHT-SCOPE SURGY MINI (NS)	THD Surgi Mini LightScope