




THD SURGY LIGHTSCOPE STERILE



Use

The SURGY LightScope anoscope is a SINGLE-USE and STERILE medical device indicated for anorectal canal analysis and surgery.

General description

The SURGY LightScope anoscope combines a practical and lightweight anoscope and a powerful integrated LED light source in one single-use device in order to give you a clear and detailed view of the inspection area.

The SURGY LightScope anoscope is ready for use and does not require any other equipment for use.

The light source is activated simply by pulling out a tab.


Features

- Made of medical-grade polycarbonate
- Complies with the biocompatibility requirements established by ISO 10993
- **Self-light, with a permanent integrated light source**
- Atraumatic insertion profile
- Extremely compact
- Ready to use
- Lightweight (40g)
- Sterilised with Ethylen Oxide
- Class I CE device
- Individually packed
- Latex-free



**Advantages:**

- Sterile device
- Does not require an additional light source
- Does not require preparation before use
- Eliminates the risk of cross-contamination
- Compact and easy to handle
- Mini-invasive and atraumatic



	Surgy mini	Surgy Maxi
E	27 mm	32 mm
LF	122 mm	122 mm
I	23	28
L	85	85

After use, dispose THD SURGY LightScope anoscopes
Like biological waste according to local laws.



*If THD SURGY LightScope anoscope need to
be disposed before it has been used (i.e. because it is expired
or damaged), do not dispose as general waste. It contains
electric components*

THD Surgi LightScope Sterile		
Code	Qty	Denomination
800124	10	THD Surgi Mini LightScope sterile
800125	10	THD Surgi Maxi LightScope sterile