




THD LIGHTSCOPE ANOSCOPE STERILE



Use

The LightScope anoscope is a DISPOSABLE and STERILE medical device deemed for investigation of the anal canal.

General Description

The LightScope anoscope can be used with the specific THD optical fiber end tip P/N 880004/880012. The LightScope anoscope has been specifically designed to be used with the THD optical fiber end tip and it has a flat white end to reflect the light and have a clear and detailed vision of the inspection field.

Available with truncated or beveled tip.

Features

- Made of a medical grade Polycarbonate
- Comply with the biocompatibility requirements of the ISO 10993 standard
- Specifically designed mandrel
- Two ends available (truncated and beveled)
- Ready to use
- Light (40 g/0.09 pound)
- CE Class Is
- Individually packaged
- Latex Free
- Sterilized by Ethylene Oxide



Advantages

- Effective lightning
- Compatible with the THD end tip (P/N. 880004/880012)
- Sterile Device

Dimensional Features	
THD LIGHTSCOPE ANOSCOPE STERILE	
Diameter (B)	25 mm (0.98 in)
Length (A)	70 mm (2.76 in)

THD LightScope Anoscope Sterile		
TRUNCATED		BEVELLED
Code	Qty	Denomination
800115 - 10	10	THD LightScope Ano Truncated Sterile
800114 - 10	10	THD LightScope Ano Bevelled Sterile