



THD ISOLATION TRANSFORMER



Use

THD Isolation Transformer is an integral component of the THD ProctoStation system.

General Description

The isolation transformer electrically isolates the main supply system from the computer which power up the peripherals.

The THD Isolation Transformer eliminates any leakage current to the connected peripherals, to comply with the required standards of patients' safety.

The THD Isolation Transformer is a certified device (IEC60601-1:2005) characterized by a very low dispersion magnetic field which does not generate disturbance and it actually protects the connected components from possible short circuits or excess voltages caused by potential malfunctions of the network where the devices are connected.



Dimensional Features	
THD ISOLATION TRANSFORMER	
Height (H)	85mm
Width (B)	150mm
Depth (T)	240mm
Weight	4.5 kg

Functional Features	
Power Absorption	300 VA
Input	115/230 V
Output	230 V on 4 outputs
Frequency	50/60 Hz
Network Switch	yes
Max. Power Output	1.3 A

ISOLATION TRANSFORMER		
Q.ty	Code	Description
1	880027	THD ISOLATION TRANSFORMER

THD360556_EN_Rev00_20180126_THD_Transformatore.docx