



THD SMOKE EVACUATOR SE02-220



Use

The Bovie Smoke Shark II is compact smoke evacuator intended to reduce the risk of electrosurgery by removing potentially harmful surgical smoke from the air.

General Description

The Bovie Smoke Shark II is a light-weight, ultra-quiet, and compact smoke evacuator that features 3 ports for maximum versatility of connections. An indicator level on the front panel allows for quickly checking your filter life. The Smoke Shark II also features a simple, push-button fan adjustment with LED indicator to easily identify the current setting. Included foot switch makes hands-free activation simple right out-of-the-box. For best results, use with Smoke-Evac Fusion Hand Piece (800206).

Technical Specifications	
Height	15.2 cm
Width	27.9 cm
Depth	39.4 cm
Weight	4.3 kgs
Line Voltage	240 VAC
Line Frequency	50-60 Hz
Noise Level	≥ 55 dBA



THD SMOKE EVACUATOR SE02-220 ACCESSORIES



THD SF35 HR REPLACEMENT FILTER P/N 800200

The 35 Hr Replacement Filter is a 4-stage filter designed to filter aerosols, gross particulate, and micro-organisms from 0.1 – 0.2 microns at an efficiency of 99.999%. A carbon filter ensures the removal and absorption of toxic gases and odors. The 3-port design accommodates a variety of tubing sizes. Item is sold separately.



THD SERS 2 REMOTE ACTIVATOR EZLINK P/N 800200

The EZLINK Remote Activator allows for automatic activation of the Smoke Shark II when connected with Electrosurgical Unit (800182) or High Frequency Desiccator (800231). Adjustable delay can be set from 0 to 10 seconds after deactivation of electrosurgical device to ensure capture of residual smoke plume. Item is sold separately.

THD SMOKE EVACUATOR SE02-220		
Q.ty	Code	Description
1	800181	THD SMOKE EVACUATOR SE02-220
1	800200	THD SF35 HR REPLACEMENT FILTER
1	800201	THD SERS 2 REMOTE ACTIVATOR EZLINK